

# Conjugates

| Code   | Quantity | Description  |
|--|----------|--|
| <b>For Monkey tissues</b>                              |          |  |
| <b>IgG FITC/Evans (M)</b>                              |          |  |
| 44847  | 3.5 mL   | Anti-(Monkey adsorbed) IgG FITC conjugate, contains Evans Blue.                |
| 44882  | 10 mL    | For use with GBMA, AICA, ATA, ASA and AACA.                                    |
| <b>IgG FITC/Evans (M) hs</b>                           |          |  |
| 44914  | 3.5 mL   | High sensitive Anti-(Monkey absorbed) IgG FITC conjugate, contains Evans Blue. |
| 44915  | 10 mL    | For use with AEA.  |
| <b>IgA FITC/Evans</b>                                  |          |  |
| 44692  | 3.5 mL   | Anti-IgA FITC conjugate, contains Evans Blue.                                  |
| 44835  | 10 mL    | For use with AEA.  |
| <b>For Rat and Mouse tissues</b>                       |          |  |
| <b>FITC/Evans (R)</b>                                  |          |  |
| 44590  | 3.5 mL   | Anti-(Rat adsorbed) IgG, IgA, IgM FITC Conjugate, contains Evans Blue.         |
| 44836  | 10 mL    | For use with ANA, AMA, ASMA, APCA and LKM.                                     |
| <b>For HEp-2, nDNA, ANCA and Rat and Mouse Tissues</b> |          |  |
| <b>IgG FITC/Evans</b>                                  |          |  |
| 44697  | 3 mL     | Anti-IgG FITC conjugate, contains Evans Blue.                                  |
| 44834  | 10 mL    | For use with HEp-2, nDNA, ANCA, ASMA, AKA, ANA, AMA, ASMA, APCA and LKM        |
| <b>Auxiliary Reagents</b>                              |          |  |
| <b>Miscellaneous</b>                                   |          |  |
| 44592  | 100 mL   | PBS (10x)  |
| 44694  | 3 mL     | Mounting medium  |
| 44846  | 10 mL    | Glycin (for GBMA)  |
| 44897  | 100 u.   | Coverslips (24 x 60 mm)  |
| 44962  | 500 mL   | PBS (10x)  |